



Recombinant Arabidopsis thaliana Membrane steroid-binding protein 2 (MSBP2)

Product Code	CSB-CF885526DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M2Z4
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MVQQIWETLKETITAYTGLSPAFFTVLALAFVYQVVSGFFVSPEVHRPRSLE VQPQSE PLPPPVLGEITEEELKLYDGSDSKPLLMAIKGQIYDVSQSRMFYGGPYAL FAGKDA SRALAKMSFEDQDLTGDISGLGAFELEALQDWEYKFMSKYVKVGTIQKKDGE GKESSEPS EAKTASAEGLSTNTGEEASAITHDETSRSTGEKIAETTEKKDVATDDDDAAKE
Source	in vitro E.coli expression system
Target Names	MSBP2
Protein Names	Recommended name: Membrane steroid-binding protein 2 Short name= AtMP2
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.