



Recombinant Rat Bcl-2-like protein 1 (Bcl2l1)

Product Code	CSB-CF002613RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53563
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSQSNRELVVDFLSYKLSQKGYWSQFSDVEENRTEAPEETEPERETPSAIN GNPSWHLA DSPA VNGATGHSSSLDAREVIPMAAVKQALREAGDEFELRYRRAFSDLTSQLH ITPGTAY QSFEQVVNELFRDGVNWGRIVAFFSFGGALCVESVDKEMQVLVSRIASWMAT YLNDHLEP WIQENGGWDTFVDLYGNNAAAESRKGQERFNRWFLTGMTVAGVVLLGSLFS RK
Source	in vitro E.coli expression system
Target Names	Bcl2l1
Protein Names	Recommended name: Bcl-2-like protein 1 Short name= Bcl2-L-1 Alternative name(s): Apoptosis regulator Bcl-X
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.