



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 37, mitochondrial (AIM37)

Product Code	CSB-CF512052SVL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C7GR78
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Sequence	MVNFYDDVDESKSHGEFPLIPVVLQNSSELSVRTIPTGNEIIESVHLTKWLRKY RNLAS QLDRYEKGWQSKIANFRLQVQHVINYSRKNIFNVDSNKHTVVPGLIALGAFF AGSIAV NRSNWGAKRLIFGHKSSILEKLCTSLPSRILLPWVLAATFKYWAPQTSQNLVN ATENDL LPADFKSYHNTWKRIYEEGYVAKKCDLKRQIDQTLQKNIRYAREQLYEKLEQ A
Source	in vitro E.coli expression system
Target Names	MIC27
Protein Names	Recommended name: Altered inheritance of mitochondria protein 37, mitochondrial
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.