



Recombinant *Saccharomyces cerevisiae* Golgi to ER traffic protein 1 (GET1)

Product Code	CSB-CF451437SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LIN5
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	MHWAAVAIAFFIVVTKFLQYTNKYHEKWISKFAPGNELSKKYLAKVKERHELKE FNNSIS AQDNYAKWTKNNRKLDSLKDEINNLIKDEIQSENKAFQAHLHLKRLLLALTPFFV FKIMYG KTPVYKLSSTSTLFPFVSGVWSQGWLYVLLHPLRTISQKWHIMEGKFGASK FDDMALQ SVSLGIWVWALMNVINGVEFIVKQLFLTPKMEAPASVETQEEKALDAVDDAIILD
Source	in vitro <i>E.coli</i> expression system
Target Names	GET1
Protein Names	Recommended name: Golgi to ER traffic protein 1 Alternative name(s): Guided entry of tail-anchored proteins 1
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.