



# Recombinant Candida dubliniensis Probable endonuclease LCL3 (LCL3)

<b>Product Code</b>	CSB-CF497031CZH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B9W9Z5
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Candida dubliniensis (strain CD36 / ATCC MYA-646 / CBS 7987 / NCPF 3949 / NRRL Y-17841) (Yeast)
<b>Sequence</b>	MPPIPPDPTEISISIFHPKVILLSAGVTTSLFFGYKFKYIRRIKTYLDLTPTIENN TK LYGYVTRVGDGDNFRFYHTPGGWIFGWGWLKIPTRKDLKDETLMIRLCGV DAPEGAHF GKPAQPFSKEALHWLREYVDGKYVTITPYSIDQYKRVARAQIWKWTGKKDIS AEMLKVG YAIVYEGKAEAEFGDNEDWYRKLESRAKLLRKGVWSLGKNLTPGEFKRIHYR GE
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	LCL3
<b>Protein Names</b>	Recommended name: Probable endonuclease LCL3 EC= 3.1.-.-
<b>Expression Region</b>	1-235
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.