



Recombinant Arabidopsis thaliana Peroxisomal membrane protein 11C (PEX11C)

Product Code	CSB-CF867987DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9LQ73
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSTLETTRAELGLVVVYLNKAEARDKICRAIQYGSKFLSDGQPGTAQNVDKNT SLARKVF RLFKFVNDLHALISVVPKGTPLPLVLLGKSKNALLSTFLFLDQIVWLGRGTGIYKD KERAE ILGRISLFCWMGSSVCTSLVEVGELGRLSASIKKLEKEIGNKDKHQNEQYRAKV EKSNER SLALIKAGMDVVVAFGLLQLAPKKVTPRVTGAFGFASLISCYQLLPSPKSKM V
Source	in vitro E.coli expression system
Target Names	PEX11C
Protein Names	Recommended name: Peroxisomal membrane protein 11C Alternative name(s): Peroxin-11C Short name= AtPEX11c
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.