



Recombinant Mouse Apoptosis regulator Bcl-2 (Bcl2)

Product Code	CSB-CF002611MO(A4)
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10417
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAQAGRTGYDNREIVMKYIHYKLSQRGYEWDAAGDADAAPLGAAPTPGIFSFQ PESNPMPA VHREMAARTSPLRPLVATAGPALSPVPPCVHLLRRAGDDFSRRYRRDFAEM SSQLHLTP FTARGRFATVVEELFRDGVNWGRIVAFFEFGGVMCVESVNREMSPLVDNIAL WMTEYLNLR HLHTWIQDNGGWDAFVELYGPSMRPLFDFSWLSLKTLLSLALVGACITLGAYL GHK
Source	in vitro E.coli expression system
Target Names	Bcl2
Protein Names	Recommended name: Apoptosis regulator Bcl-2
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.