



Recombinant *Oryza sativa* subsp. japonica Peroxisomal membrane protein 11-1 (PEX11-1)

Product Code	CSB-CF604686OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10SM7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MSTLDATRAELGLVVLNKA EARDKICRAIQY GSKFISNGQP GTAQDVDRSTT LARKVF RLLKWN DLHGLISPPAKGTPL TLVLLGKSKNALLSTFLFDQFVWLGRTGIYK NKERTD RIVRISLYCWMASSVCAGLVELGELKRLSKSMRKLARELRD TDKYENDQYKSK MKQSDER LLALVKAAMDVVAVGLLQLSPK KITPRVTGAFGFV TSLISCYQQLPSRAPAIKV KA
Source	in vitro <i>E.coli</i> expression system
Target Names	PEX11-1
Protein Names	Recommended name: Peroxisomal membrane protein 11-1 Alternative name(s): OsPEX11-1 Peroxin-11-1
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.