



Recombinant Arabidopsis thaliana Ubiquitin-conjugating enzyme E2 34 (UBC34)

Product Code	CSB-CF874417DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SHI7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MAEKACIKRLQKEYRALCKEPVSHVVARPSPNDILEWHYVLEGSEGTPFAGGF YYGKIKF PPEYPYKPPGITMTTPNGRFMTQKKICLSMSDFHPESWNPMSVSSILTGLLS FMMDTSP TTGSVNTTVIEKQRLAKSSLAFNCKTPAFRKLFPYVEKYNQQLAEQATTQLT TPESPQ KSDTKVESEKTIDPTKGDSEGLKERKKNKQGLPAWIILLVSVFGVVMALPL LQL
Source	in vitro E.coli expression system
Target Names	UBC34
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 34 EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein 34
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.