



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL57 (ATL57)

Product Code	CSB-CF890362DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SJJ7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	METKHLRKLLQLYQACGGEQELTKTVQNVTSPLTPPPPQPPSALDSTMALT IFILLVA LFFMGFFSVYFRHFADSSSSTVDISSMPRTRSSRMSPRRLSTSVVVSRPYSFR RGLDSQA VRSLPVYRYTKAAKQRNEDCVICLSDFEEGETVKVIPHCGHVFHVDCVDTWLS SYVTCPL CRSNQLFSDKDLGMQEPPDQDSAEHDTCDGVDTCVRRCSSCSLGLQRTGL ERSLSL
Source	in vitro E.coli expression system
Target Names	ATL57
Protein Names	Recommended name: RING-H2 finger protein ATL57
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.