



Recombinant Zea mays Cell number regulator 6 (CNR6)

Product Code	CSB-CF467294ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6SGC5
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MAEDATSSHPSRYVKLTKDQDAPAEDIRPGELNQPVHVPQLEGRRRCSECGQV LPESYEPP ADEPWTTGIFGCTDDPETCRTGLFCPCVLFGRNVEAVREDIPWTTPCVCHAVF VEGGITL AILTAIFHGVDPRTSFLIGEGLVFSWWLCATYTGIFRQGLQRKYHLKNSPCDPC MVHCCCL HWCANCQEHRERTGRLAENNAVPMTVVNPPPQEMSMLEEVEEKGAEKSEH DDVEVIPL
Source	in vitro E.coli expression system
Target Names	CNR6
Protein Names	Recommended name: Cell number regulator 6 Alternative name(s): ZmCNR06
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.