



Recombinant Human Apoptosis regulator Bcl-2 (BCL2)

Product Code	CSB-CF002611HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10415
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAHAGRTGYDNREIVMKYIHYKLSQRGYEWDA GDVGAAPPGAAPAGIFSSQ PGHTPHPA ASRDPVARTSPLQTPAAPGAAAGPALSPVPPVVHLTLRQAGDDFSRRYRRDF AEMSSQLH LTPFTARGRFATVVEELFRDGVNWGRIVAFFEFGGVMCVESVNREMSPLVDNI ALWMTEY LNRHLHTWIQDNGGWDAFVELYGPSMRPLDFDFSWLSLKTLLSLALVGACITLG AYLGHK
Source	in vitro E.coli expression system
Target Names	BCL2
Protein Names	Recommended name: Apoptosis regulator Bcl-2
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.