



Recombinant Rhizobium meliloti Protein exoD (exoD)

Product Code	CSB-CF700601RKU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q52923
Product Type	Transmembrane Protein
Immunogen Species	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Sequence	MVCQARFGVKNASRGHRMERPQTVKGGKIMAVEFGDSQSRSLSDTLTGMIASIR GNTITLRE LMIEIGEQQGFLLLCALLTLPFLIPVSIPGVSTVFGAAIILISLAILNRMPWLPKRILD R EIAATEKLVPTLRKGAALVSKLDRYVRPRLNFLTEGALMNRFNGLMIMAGGVLLM FPLGLI PLSNTLPGIAILLLSLGGIQRDGLMVAGGYFFLVATTVYFAVLGYAAFAAGQGLS HFFVS
Source	in vitro E.coli expression system
Target Names	exoD
Protein Names	Recommended name: Protein exoD
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.