



Recombinant Arabidopsis thaliana Uncharacterized ribosomal S3-like protein AtMg00690, mitochondrial (AtMg00690)

Product Code	CSB-CF308992DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92511
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MRSSVLRSLRGLVINLESTRKLRSLRTNIVPGRKKGQKSIKSKNMARKGNPIL VRLGKN RSSDSSWFSAEALLGCLYFFIYFVAPTLGPVLFLLRLIHFWGLRLGLGNENFH FGVGPD GGATGLDLNQPPQEQQPTLGVNRAALDLNELPPVHLLYAEVEGPQSTKAQND VMLAHLNQ VQNLTRDLQTEPNIWRRQALIDILDWEVRSLQRHFRIFRQRDRLREVQRSWLR EQLNRYR
Source	in vitro E.coli expression system
Target Names	AtMg00690
Protein Names	Recommended name: Uncharacterized ribosomal S3-like protein AtMg00690, mitochondrial Alternative name(s): ORF240a
Expression Region	1-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.