



Recombinant Human B-cell receptor-associated protein 29 (BCAP29)

Product Code	CSB-CF002591HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UHQ4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTLQWAAVATFLYAEIGLILIFCLPFIPPQRWQKIFSFNVWGKIATFWNKAFLLI LLI VLFLDAVREVRKYSSVHTIEKSSTSRPDAYEHTQMKLFRSQRNLYISGFSLFFW LVLRRRL VTLITQLAKELSNKGVLKTQAENTNKAACKFMEENEKLRILKSHGKDEECVLE AENKKL VEDQEKLKTELKRTSDALSKAQNDVMEMKMQSERLSKEYDQLLKEHSELQDR LERGNKKR L
Source	in vitro E.coli expression system
Target Names	BCAP29
Protein Names	Recommended name: B-cell receptor-associated protein 29 Short name= BCR-associated protein 29 Short name= Bap29
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.