



Recombinant Mouse Peroxisomal membrane protein 11C (Pex11g)

Product Code	CSB-CF744064MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6P6M5
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MALLNRLASALESHRVRDRLIRTLGYCCQLIGGVLVEQCPNRSEVGRLLVVS AQFNHCR TVLRLFDDLAMFVYTKQYGLGTKEEDIFIRWLSVLSNVTDQLYYPCEHIAWAAD AKVLRV DSAWWWTLNLTALWTLSSLLGAVKALWTMLKLRQKLRSPGTGTSASQLPRSKRR AMEARICS EVLTLNLSNLADLANAVHWLPRGVLWAGRFPPWLVGLMGTISSILSTCQAVRAG RQAEADS P
Source	in vitro E.coli expression system
Target Names	Pex11g
Protein Names	Recommended name: Peroxisomal membrane protein 11C Alternative name(s): Peroxin-11C Peroxisomal biogenesis factor 11C Protein PEX11 homolog gamma Short name= PEX11-gamma
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.