



Recombinant Human cytomegalovirus Envelope glycoprotein UL132 (UL132)

Product Code	CSB-CF304806HWV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69338
Product Type	Transmembrane Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Sequence	TNMTSSTNVPTSTSSRNTVESTTSSEPTTETNMTTARESSVHDARNDEIMKVL AILFYIV TGTSIFSFI AVLIAVVYSSCCKHPGRFRFADEEAVNLLDDTDDSGGSSPFGSGS RRGSQI PAGFCSSSPYQRLETRDWDEEEEEASAARERMKHDPENVIYFRKDGNLDTSFV NPNYGRGS PLTIESHLSDNEEDPIRYYSVYDELTA SEMEEPSNSTSWQIPKLMKVAMQPV SLRDPEY D
Source	in vitro E.coli expression system
Target Names	UL132
Protein Names	Recommended name: Envelope glycoprotein UL132 Short name= L3
Expression Region	30-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.