



Recombinant *Oryza sativa* subsp. *indica* Peroxisomal membrane protein 11-3 (PEX11-3)

Product Code	CSB-CF386076OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2XFQ8
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MAAAAAAAGSSDSRKPAAHPPPRDFLVHVEAYLSRRDGVDKLLKISRYAARLA LAAGPLP PAASARLKSFEVVGLSRKAFRLGKVFQNVNALRAHPHPPPAVALLAYGGEG VYYFLEQF VWLAKAGLLPAHLLPRLQRLSAWAELLYVGSITIKLEEIGKLESSVKMRLKEG CREESD VVRTLRVKLLKRMSSVVQDVADAVMALGDVTDGKGLLGSSTLMASAGLLSALI SAHKNWN SC
Source	in vitro <i>E.coli</i> expression system
Target Names	PEX11-3
Protein Names	Recommended name: Peroxisomal membrane protein 11-3 Alternative name(s): OsPEX11-3 Peroxin-11-3
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.