



Recombinant Danio rerio Uncharacterized protein C18orf19 homolog A (si:ch211-105d11.2)

Product Code	CSB-CF629887DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1MTD4
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MQRFVCSGVFRQTSVLLQISPVSHCHTRRMQDSPSHRSLKQLISSQRWLNST TALCSPTS QKTSASGQEEPDPPLHDKSSGLIQRFKKTFKQYGKVMIPVHLLTSTMWFGTFYY AAMKGVN LVPFLEYVGFDPKVVKLENSQSGYALTAYAMYKIATPARYTVTLGGTSLSVKY LRKHGY MSTPPPVKEYLQEKMEETKERISGKMEETKDRFSERMEETKDKFNEKLQETK DKVSFRKK KE
Source	in vitro E.coli expression system
Target Names	si:ch211-105d11.2
Protein Names	Recommended name: Uncharacterized protein C18orf19 homolog A
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.