



Recombinant Arabidopsis thaliana E3 ubiquitin-protein ligase RMA3 (RMA3)

Product Code	CSB-CF815421DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8GUK7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MEGNFFIRSDAQRADNGFIAKQKPNLTTAPTAGQANESGCFDCNICLDTAHD PVVTLCGHLFCWPCIKWLHVQLSSVSVDQHQNNCPVCKSNITITSLVPLYGR GMSSPSSTFGSKKQDALSTDIPRRPAPSALRNPITSASSLNPSLQHQTLSPSFH NHQYSPRGFTTTESTDLANAVMMSFLYPVIGMFGDLVYTRIFGFTNTIAQPYQ SQRMMQREKSLNRVSIFFLCCILCLLLF
Source	in vitro E.coli expression system
Target Names	RMA3
Protein Names	Recommended name: E3 ubiquitin-protein ligase RMA3 EC= 6.3.2.- Alternative name(s): Protein RING membrane-anchor 3
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.