



Recombinant Biopolymer transport protein exbB (exbB)

Product Code	CSB-CF359669EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ABU8
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MGNNLMQTDLSVWGMYPHADIVVKCVMIGLILASVVTWAIFFSKSVEFFNQKR RLKREQQ LLAEARSLNQANDIAADFGSKSLSLHLLNEAQNELELSEGSDDNEGIKERTSFR LERRVA AVGRQMGRGNGYLATIGAI SPFVGLFGTVWGIMNSFIGIAQTQTTNLAVVAPGI AEALLA TAIGLVAAIPAVVIYNV FARQIGGFKAM LGDVAAQVLLLQSRDL DLEASAAAHPV RVAQK LRA G
Source	in vitro E.coli expression system
Target Names	exbB
Protein Names	Recommended name: Biopolymer transport protein exbB
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.