



Recombinant Arabidopsis thaliana Reticulon-like protein B7 (RTNLB7)

Product Code	CSB-CF868101DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M145
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MEEEKLEIVGPLEEPLMGNIVPEEINGLDSLTSDDSDSEKPDSPVPINAPIYRMF GRERP IHMVLGGGKPADVLLWRDKKVTLGLLSAVTVIWLFLFGFGGRLLTSLCRGSILF LLLSFL WSNALNKSPENMMDIYIPEKPLLQAASAMTFELNCAFATLRSIALERDIKNFVM AVIGLW LVSVIGNWFSFLSLLYICFVLIHTVPMLYEKYEDEIDPIAEKAVIEMKKHYQVFEA KFLS KIPH
Source	in vitro E.coli expression system
Target Names	RTNLB7
Protein Names	Recommended name: Reticulon-like protein B7 Short name= AtRTNLB7
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.