



Recombinant Human B-cell receptor-associated protein 31 (BCAP31)

Product Code	CSB-CF002592HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51572
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SLQWTAVATFLYAIEVFVLLLCIPFISPKRWQKIFKSRLVELLVSYGNTFFVVLIV ILVL LVIDAVREIRKYDDVTEKVNQLQNNPGAMEHFHMKLFRAQRNLYIAGFSLLSFL LRRLVT LISQQATLLASNEAFKKQAESASEAAKKYMEENDQLKKGAAVDGGKLDVGNA EVKLEEN RSLKADLQKLKDELASTKQKLEKAENQVLAMRKQSEGLTKEYDRLLLEEHAQLQ AAVDGPM DKKEE
Source	in vitro E.coli expression system
Target Names	BCAP31
Protein Names	Recommended name: B-cell receptor-associated protein 31 Short name= BCR-associated protein 31 Short name= Bap31 Alternative name(s): 6C6-AG tumor-associated antigen Protein CDM p28
Expression Region	2-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.