



Recombinant Dictyostelium discoideum Cytochrome c1, heme protein, mitochondrial (cyc1)

Product Code	CSB-CF706679DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54D07
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	ITTSTLLSDDFVQPTSYPWEHRFPWQSYDHAAIRRGHQVYTQVCSTCHSLNL ISYRHLA NVAYTPEELKAMAADTTVMDGPDSEGDMFERKGGQVTDHFHPKPYPNSNAARF ANNALPPD LSLVIKARGAHEDYVFSLLTGYCEPPAGVRVVGQYFNPYFPGTKIAMAPPLA DGMVEYD DGTDNSMSQMAKDVSTFLCWASEPEHDDRKCLGMKVLLGVSIIALPLFYWKR LKWSVIKT RKISFKD
Source	in vitro E.coli expression system
Target Names	cyc1
Protein Names	Recommended name: Cytochrome c1, heme protein, mitochondrial Short name= Cytochrome c-1 Alternative name(s): Complex III subunit 4 Complex III subunit IV Cytochrome b-c1 complex subunit 4
Expression Region	29-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.