



Recombinant Human Peroxisomal membrane protein 11A (PEX11A)

Product Code	CSB-CF017795HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75192
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDAFTRFTNQQTQGRDRLFRATQYTCMLLRYLLEPKAGKEKVVMMKLLKLESSV STGRKWFR LGNVVHAIQATEQSIHATDLVPRCLCLLANLNRVIYFICDTILWVRSVGLTSGINK EKWR TRAAHHYYYSLLLSLVRDLYEISLQMKRVTCDRAKKEKSASQDPLWFSVAEEE TEWLQSF LLLLFRSLKQHPPLLLDTVKNLCDILNPLDQLGIYKSNPGIIGLGLVSSIAGMITV AYP QMCLKTR
Source	in vitro E.coli expression system
Target Names	PEX11A
Protein Names	Recommended name: Peroxisomal membrane protein 11A Short name= HsPEX11p Alternative name(s): 28 kDa peroxisomal integral membrane protein Short name= PMP28 Peroxin-11A Peroxisomal biogenesis factor 11A Protein PEX11 homol
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.