



Recombinant Human CKLF-like MARVEL transmembrane domain-containing protein 2 (CMTM2)

Product Code	CSB-CF844708HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TAZ6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPKAAKGAKPEPAPAPPPPGAKPEEDKKDGKEPSDKPQKAVQDHKEPSDK PQKAVQPKH EVGTRRGCRRYRWELKDSNKEFWLLGHAEIKIRSLGCLIAAMILLSSLTVHPILR LIITM EISFFSFFILLYSFAIHRYIPFILWPISDLFNDLIACAFLVGAVVFAVRSRRSMNLH YLL AVILIGAAGVFAFIDVCLQRNHFRGKKAKKHMLVPPPGKEKGPQQGKGPEPAK PPEPGKP PGPAGKGG
Source	in vitro E.coli expression system
Target Names	CMTM2
Protein Names	Recommended name: CKLF-like MARVEL transmembrane domain-containing protein 2 Alternative name(s): Chemokine-like factor superfamily member 2
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.