



Recombinant Pongo abelii Cytochrome c oxidase assembly protein COX11, mitochondrial (COX11)

Product Code	CSB-CF728956PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R7U6
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	AERVEPFLRPECSGTGGAGRGLRWLGTWKRCSLGARHPALQPPRRHKNSNP FTRAQKEEW RRRNKTTLTYYAAVAVGMLGASYAAVPLYRLYCQTTGLGGSVAVAGHASDKIEN MVPVKDR IIKISFNADVHASLQWNFRPQQTEIYVVPGETALAFYRAKNPTDKPVGISTYNIV PFEA GQYFNKIQQCFCFEEQRLNPQEEVDMPVFFYIDPEFAEDPRMIKVDLITLSYTF FEAKEG HKLPVPGYN
Source	in vitro E.coli expression system
Target Names	COX11
Protein Names	Recommended name: Cytochrome c oxidase assembly protein COX11, mitochondrial
Expression Region	29-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.