



Recombinant Vaccinia virus Protein L1 (L1R)

Product Code	CSB-CF324949VAA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20540
Product Type	Transmembrane Protein
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Sequence	GAAASIQTTVNTLSERISSKLEQEANASAQTKCDIEIGNFYIRQNHGCNLTVKN MCSADA DAQLDAVLSAATETYSGLTPEQKAYVPAMFTAALNIQTSVNTVVRDFENYVKQ TCNSSAV VDNKLKIQNVIIDECYGAPGSPTNLEFINTGSSKGNCAIKALMQLTTKATTQIAP RQVAG TGVQFYMIVIGVILAALFMYYAKRMLFTSTNDKIKLILANKENVHWTTYMDTFFR TSPM VIATTDMQN
Source	in vitro E.coli expression system
Target Names	L1R
Protein Names	Recommended name: Protein L1 Alternative name(s): Virion membrane protein M25
Expression Region	2-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.