



Recombinant Schizosaccharomyces pombe UPF0494 membrane protein PB2B2.07c (SPBPB2B2.07c)

Product Code	CSB-CF862080SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9HDU8
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MSNPESAKKQVDPPGYNELFMVRDTRNVDLRLERGLELCKPEKVNKQNLFTNIIK PQKDKIN IKTDKIKFFLNLFTEFSKFHDSCYPDGRISTRSLRWPLLIWLCILIVFAIDKNFE VKD FLSIWINESFINENRFYSEIWGPPIAYICLFILLLLGLICMFPLHLRCRVCVLALRETG MI IAVLGAALGMIIAALGATITGLLYFSHWALYKVVILALDLKIEPFKDEIAFLTLPTHN GE TLSIRDNNQS
Source	in vitro E.coli expression system
Target Names	SPBPB2B2.07c
Protein Names	Recommended name: UPF0494 membrane protein PB2B2.07c
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.