



Recombinant Human Receptor expression-enhancing protein 2 (REEP2)

Product Code	CSB-CF887109HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BRK0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MVSWIISRLVVLIFGTLYPAYSSYKAVKTKNVKEYVKWMMYWIVFAFFTTAETL TDIVLS WFPFYFELKIAFVIWLLSPYTKGSSVLYRKVFVHPTLSNKEKEIDEYITQARDKSY ETMMR VGKRGLNLAANAAVTAAAKGVLSEKLRFSMQDLTLIRDEDALPLQRPDGRRLR PSPGSL DTIEDLGDDPALSLSSTNPADSRTEASEDDMGDKAPKRAKPIKKAPKAEPLA SKTLKTR PKKKTSGGGDSA
Source	in vitro E.coli expression system
Target Names	REEP2
Protein Names	Recommended name: Receptor expression-enhancing protein 2
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.