



Recombinant Human Uncharacterized protein C1orf43 (C1orf43)

Product Code	CSB-CF874824HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BWL3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASGSNWLSGVNVVLMAYGSLVFVLLFIFVKRQIMRFAMKSRRGPHVPVGH NAPKDLKE EIDIRLSRVQDIKYEPQLLADDDARLLQLETQGNQSCYNYLYRMKALDAIRTSEI PFHSE GRHPRSLMGKNFRSYLLDLRNTSTPFKGVRKALIDTLLDGYETARYGTGVFGQ NEYLRYQ EALSELATAVKARIGSSQRHHQSAAKDLTQSPEVSPTTIQVTYLPSSQKSKRAK HFLELK SFKDNyntLESTL
Source	in vitro E.coli expression system
Target Names	C1orf43
Protein Names	Recommended name: Uncharacterized protein C1orf43 Alternative name(s): Hepatitis C virus NS5A-transactivated protein 4 Short name= HCV NS5A- transactivated protein 4 Protein NICE-3 S863-3
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.