



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL28 (ATL28)

Product Code	CSB-CF740632DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6NKR1
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSSTTSIDSIPATAVFPSVSMPVTVVLTVGVLVIFAGFFSLFLWQFLLNRLFTT WNLQR TPYGDLIHVATPPENTGLDPFIIRSFPVFHYSSATKKNHGTECAICLSEFSDEDT VRLIT VCRHPFHNSCIDLWFELHKTCPVCRCELDPGMIGSGRLESFHNTVTITIQDINH DEENPP TAGSSKRLIEASAWRFSRSHSTGHFMVKTTDANVKSRRHYQTGSCVSFDEL TRYEGAGW QWLGDSSHISRIEV
Source	in vitro E.coli expression system
Target Names	ATL28
Protein Names	Recommended name: RING-H2 finger protein ATL28
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.