



Recombinant Human Emerin (EMD)

Product Code	CSB-CF007633HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50402
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDNYADLSDELTTLLRRYNIPHGPVVGSTRRLYEKKIFEYETQRRRLSPPSSS AASSYS FSDLNSTRGDADMYDLPKKEDALLYQSKGYNDYEEESYFTTRTYGEPESAG PSRAVRQS VTSFPDADAFHHQVHDDDLLSSSEEECKDRERPMYGRDSAYQSITHYRPVSA SRSSLDLS YYPTSSSTSFMSSSSSSSSWLTRRAIRPENRAPGAGLGQDRQVPLWGQLLLF LVFVIVLF FIYHFMQAEEGNPF
Source	in vitro E.coli expression system
Target Names	EMD
Protein Names	Recommended name: Emerin
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.