



Recombinant Arabidopsis thaliana Reticulon-like protein B5 (RTNLB5)

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF525760DOA |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O82352 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Sequence | MAEEIEKSVPTTEESLMEKISEKIHHDSSSSSESEYEKPDSPSAVKAKIYRMFG REKPVH KVLGGGKPADVFLWRDKKLSGAVLGVATAIWVLFELVEYHLLSLLCHISILALG GLFLWS NAHTLINKTSPQIPEIHVPEEAFLVVASSLRNELNQAFVILRSIALGRDLKKFLMV VVGL WIISVVGWVFNFLTLVYICFVILHTVPMLEYKHEDKVDPLAEKAMKELQKQYVV FDEKVL SKIPIASLKAKAKLG |
| Source | in vitro E.coli expression system |
| Target Names | RTNLB5 |
| Protein Names | Recommended name: Reticulon-like protein B5 Short name= AtRTNLB5 |
| Expression Region | 1-255 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |