



Recombinant Arabidopsis thaliana Reticulon-like protein B4 (RTNLB4)

Product Code	CSB-CF884223DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FFS0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MVEDHKHEESILEKIVEKIHGHGDSSSLSDSDDDKKSTSSSSSSFKSKIYRLFG REKPVH KVLGGGKPADIFLWRNKKVSGGVLGAVTASWVLFELFEYHLLAFLCHFAIFALA ALFLWS NACTFIHKSTPHIPEVHIPEDPILQLVSGLRIEINRGLTLRNIASGKDVKKFILVIA GL WVLSIIGSCYNFLTIFYTATVLLFTIPVLYEKYEDKVDAYGEKAMREIKKQYAVL DEKVL RKVISKIPRGALNKKKD
Source	in vitro E.coli expression system
Target Names	RTNLB4
Protein Names	Recommended name: Reticulon-like protein B4 Short name= AtRTNLB4 Alternative name(s): VirB2-interacting protein 3
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.