



Recombinant Pongo abelii Peroxisomal membrane protein 11B (PEX11B)

Product Code	CSB-CF732819PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5RF10
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	MDAWVRFSAQSQARERLCRAAQYACSLLGHVLRHGASPELQKQIRQLESHL SLGRKLLR LGNSADALESAKRAVHLSDVVLRFCITVSHLNRALYFACDNVLWAGKSGLAPR VDQEKWA QRSFRYYLFLIMNLSRDAYEIRLLMEQESSACSRRLKGGGGVPGGSETGGL GGPGTPG GHLPQLALKLRLQVLLARVLRGHPPLLLDVVRNACDLFIPLDKLGLWRCGPGI VGLCGL VSSILSILTLIYPWLRLKP
Source	in vitro E.coli expression system
Target Names	PEX11B
Protein Names	Recommended name: Peroxisomal membrane protein 11B Alternative name(s): Peroxin-11B Peroxisomal biogenesis factor 11B Protein PEX11 homolog beta Short name= PEX11-beta
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.