



Recombinant *Drosophila melanogaster* RPII140-upstream gene protein (140up)

Product Code	CSB-CF305633DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81928
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MNFLWKGRRFLIAGILPTFEGAADIVDKENKTYKAFLASKPPEETGLERLKQM FTIDEF GSISSELNSVYQAGFLGFLIGAIYGGVTQSRVAYMNFMENNQATAFKSHFDAK KKLQDQF TVNFAKGGFKWGWVRVGLFTTSYFGIITCMSVYRGKSSIYEYLAAGSITGSLYKV SLGLRG MAAGGIIGGFLGGVAGVTSLLLMKASGTSMEEVRYWQYKWRLDRDENIQQAF KKLTEDEN PELFKAHDEKTSEHVSLDTIK
Source	in vitro <i>E.coli</i> expression system
Target Names	140up
Protein Names	Recommended name: RPII140-upstream gene protein
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.