



Recombinant *Saccharomyces cerevisiae* Genetic interactor of prohibitin 7, mitochondrial (GEP7)

Product Code	CSB-CF457258SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LHC6
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	AVAGNLLVKRFYQPKLERIPPASLLLKQKIRLAQNGSTTSTENPISFSQTMSEIF SVLQP SAPDLDEDETSGLKRDHLLTERLNNGELGVIMNKFFNPSSTHNNQLIDTNILLQ NFPKLS GNDLDDLLDFAIN EKMRGNWNDLKQDFIQLWYYKSFGLGPRTQFVL TNSSPSL RSQFLKL PFIEYNWFL LQNNKNANILP ADVQNVVKV FHLDDKRF SWKSIDPF SKAIISFVVF VSIYV WLDESAKQKTKELPAQKSTVISE
Source	in vitro <i>E.coli</i> expression system
Target Names	GEP7
Protein Names	Recommended name: Genetic interactor of prohibitin 7, mitochondrial
Expression Region	25-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.