



Recombinant Escherichia coli Uncharacterized protein yjfZ (yjfZ)

Product Code	CSB-CF340081ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P39308
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MTLPTTIYSFPAYLSRFSSTDKPVKLFHQYARATLLSNRGRDHNC DGRRTVEI HKLDLS DWQAFNKLATRCNAYDGITMNGDNSFGWNHEATLDNIHAQKYNKAYAGARLT AELKYLLQ DVESFEPNSKYTIHEVVLGPGYGTPDYTGQTIGYVVTLPAQMPNCWSSSELPTI DLYIDQL RTVTGVSNALGFIIAALLNAYS DLP HDL KIGLRSLSSSAAIYSGLGFERVPQERDI SCAR MYLTPANHPDLWTQENGEWIYLRN
Source	in vitro E.coli expression system
Target Names	yjfZ
Protein Names	Recommended name: Uncharacterized protein yjfZ
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.