



Recombinant Danio rerio Receptor expression-enhancing protein 2 (reep2)

Product Code	CSB-CF679878DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4KMI4
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MVSWIISRMVVLAFGTLYPAYSSYKAVKTKNVKEYVKWMMYWIVFALFTTAETI TDMLLS WFPFYFELKIAFVIWLLSPYTKGSSVLYRKFVHPTLSNKEREIDEYITQAKDRSY DTMMR FGRRGLNIAATAAVTAATKGQGVLSSEKLRFSMQDLTLIQNEDELQLEGGDDT HTAATLP RAKTATRTVRATPVPADTESQHSSRSDDQSDSRTEHSDEDAADKAPKRIAIR AAKKPAA AKTEQTTKTVKKAPKKKPTTANNVAESP
Source	in vitro E.coli expression system
Target Names	reep2
Protein Names	Recommended name: Receptor expression-enhancing protein 2
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.