



Recombinant Mouse Integral membrane protein 2C (Itm2c)

Product Code	CSB-CF011906MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q91VK4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MVKISFQPAVAGIKADKADKAAASGPASASAPAAEILLTPAREERPPRHRSRKG GSVGGV CYLSMGMVVLLMGLVFASVYIYRYFFLAQLARDNFFHCGVLYEDSLSSQIRTRL ELEEDV KIYLEENYERINVPVPQFGGGDPADIIHDFQRGLTAYHDISLDKCYVIELNTTIVL PPRN FWELLMNVKRGTYLPQTYIIQEEMVVTEHVDRKEALGSFIYHLCNGKDTYRLR RRSTRRR INKRGGKNCNAIRHFENTFVVETLICGVV
Source	in vitro E.coli expression system
Target Names	Itm2c
Protein Names	Recommended name: Integral membrane protein 2C Alternative name(s): Transmembrane protein BRI3 Cleaved into the following chain: 1. CT-BRI3
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.