



Recombinant Human Protein ARV1 (ARV1)

Product Code	CSB-CF867134HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H2C2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGNGGRSGLQQGKGNVDGVAATPTAASASCQYRCIECNQEAKELYRDYNHG VLKITICKS CQKPVDKYIEYDPVILINAILCKAQAYRHILFNTQINIHGKLCIFCLLCEAYLRWW QLQ DSNQNTAPDDLIRYAKEWDFYRMFAIAALEQTAYFIGIFTFLWVERPMTAKKKP NFILL KALLSSYGKLLLIPAVIWEHDYTSVCLKLIKVFVLTSNFQAIRVTLNINRKLSFLA VLS GLLLESIMVYFFQSMEDVVGSDYAIFKSQDF
Source	in vitro E.coli expression system
Target Names	ARV1
Protein Names	Recommended name: Protein ARV1 Short name= hARV1
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.