



Recombinant Human Probable E3 ubiquitin-protein ligase RNF217 (RNF217)

Product Code	CSB-CF844717HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TC41
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKVQLGQVEIKCPITECFEFLEETTVVYNLTHEDSIKYKYFLELGRIDSSTKPCP QCKHF TTFKKKGHIPTPSRSESKYKIQCTCQFVWC FKCHSPWHEGVNCKEYKKGDK LLRHWASE IEHGQRNAQKCPKCKIHIQRTEGCDHMTCSQCNTNFCYRCGERYRQLRFFGD HTSNLSIF GCKYRYLPERPHLRRLVRGSVCAGKLFIAPLIMVLGLALGAIWVVEEIKTYWN LISGR TRNQTQHLAPQPVLLSDMLYCLKQVFICISYLLPL
Source	in vitro E.coli expression system
Target Names	RNF217
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase RNF217 EC= 6.3.2.- Alternative name(s): IBR domain-containing protein 1 RING finger protein 217
Expression Region	1-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.