



Recombinant Human Uncharacterized protein C17orf78 (C17orf78)

Product Code	CSB-CF847630HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N4C9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDTILVFSLIIASYDANKKDLRDSSCRLEQLPGIFPKDVRSIRELQMQETHETK RTTFI QNRTIATLQCLGSDSKVKVNLVYLERRPKVKHILKNLRIIAAPRRNSSASSSCHL IPTSK FQTGSELLKGKAFLPGISQCKVLGASSETFPTTAPSITPGNKEGEKTTSTDTDEN LEKRQK WSIVVKILIAVTLTLLSGVAIVFVIFEVPCPYQCLGARKLCQCQWLWRWQKKGG QPPGTA ESKPDSQPQKVGQDAANSSNPKAAEITVIHQTYF
Source	in vitro E.coli expression system
Target Names	C17orf78
Protein Names	Recommended name: Uncharacterized protein C17orf78
Expression Region	1-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.