



Recombinant Human Transmembrane protein 53 (TMEM53)

Product Code	CSB-CF023853HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6P2H8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASAELDYTIEIPDQPCWSQKNSPSPGGKEAETRQPVVILLGWGGCKDKNLA KYS AIYHK RGCIVIRYTAPWHMVFFSESLGIPSLRVLAQKLELLFDYEIEKEPLL FHVFSNG GVMLY RYVLELLQTRRFCRLRVVGTIFDSAPGDSNLV GALRALAA ILERRAAMLRLLLL V AFALV VVL FHVLLAPITALFH THFYDRLQDAGSRWPELYLYSRADEVVLARDIERMVEA RLARRV LARSVDFVSSAHVSHLRDYPTYTSLCVDFMRNCVRC
Source	in vitro E.coli expression system
Target Names	TMEM53
Protein Names	Recommended name: Transmembrane protein 53
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.