



Recombinant Danio rerio N-acetyltransferase 14 (nat14)

Product Code	CSB-CF603901DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0P4A4
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MVRLDLADVVLRRMQEKDIEAVKALIKEGCEGTENRLLHLLTRPLALLLLAILSS ILRC VLHSFVLALVIPVFISVIYLKLTIPRSAGILGSCRPYWDYIGSSYHADTEPDLPNP HLGR AKLTTNQEKTRRRKKAKEKEKMNESQVDEDELKQRAKVAGEVWVADSDGEI VGCVARDG WSRDGVCRVCRLLVVQCWYRREGLGRLLVQGLESRTKQKGVCRVYAHVPIPS KVGEAFFRR LGYRLQGETAGIEEEEEEDDYEDPEKGWLGYP LTKVFKDL
Source	in vitro E.coli expression system
Target Names	nat14
Protein Names	Recommended name: N-acetyltransferase 14 EC= 2.3.1.-
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.