



Recombinant Danio rerio Uncharacterized protein C18orf19 homolog B (zgc:113036)

Product Code	CSB-CF707702DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5CZQ0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<pre>MQRLWAPVTLRRVLLLRSVYLPHSWAMQEPRLAVLSPRYFSCNRAIRAKEAH KTSTEEQE EVPLNPPQSLAGTEGLYKADSEPVPHNKGIDIDPLQDKSIGIFQRFKTKFKQYGK VMVPVH IVTSTVWFGSFYYAAMKGVNLPFLEFIGLPDWIVGILRDSQGGYALTAYAMYK LATPAR YVTMGGTSLSVQYLRKHGYLSTPPPVKFLQDKMEETRELLTEKMEETKERF SEKMEET KELLSERMEETKERFSEKDKFSEKDKLQETKDKMSFRKKAD</pre>
Source	in vitro E.coli expression system
Target Names	zgc:113036
Protein Names	Recommended name: Uncharacterized protein C18orf19 homolog B
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.