



Recombinant Danio rerio Cholesterol 25-hydroxylase-like protein 1, member 1 (ch25hl1.1)

Product Code	CSB-CF710054DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q567X1
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MWNISEVVFQLPTSSASDRVLQPLWDYLLLRHYTLISSPFFPVLLAFSSYIIFSV PFAVL DVLGEKAPLFKYKIQKDRRPTVGMMLRTLWTAVYNHLVFLPAVLITNVVMPM PPLPTVA PTVWEMFSGGLGALLVFDTQYFLWHMVHHKNPHLYRWWHAIHHDYISPFWS TQHLSGVE LMTVGFWSNIDPILLKCHPLTVWTLTVYSIWMSVEDHIGYDLPFSPGHLVPFGL LGGAMA HDMHHQKPSSNFAPFFSHWDKIFGTAITVKLTQKSEKEKQVA
Source	in vitro E.coli expression system
Target Names	ch25hl1.1
Protein Names	Recommended name: Cholesterol 25-hydroxylase-like protein 1, member 1 EC=1.14.99.-
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.