



Recombinant *Saccharomyces cerevisiae* Pheromone-regulated membrane protein 4 (PRM4)

Product Code	CSB-CF622487SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q12498
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MIADSSVLKKHTAIKRSTRIISLTLVLLGVFSLLLLTWNDSLEFYNSADPSENKKN SEEE SEKKFVYKLPNLLKTADSFLSNENELNFQVKKEEISNIQSEVEVDIPEPSSKATS KFSSR SFQTDNVVTATTTTTLNPRSSSLALQKNCDHKKFDPRTDFLDIIRTSPAFLIKS SQADS IFLKNLLQREFEISPELATVDLEKHSYGYELEKYIKQNKLNIDTSAALESIQSPYL FLNG ISVINRGMVRDIIEPHSGLLLLPLLKSEARGNLLVEKKDIPSNS
Source	in vitro <i>E.coli</i> expression system
Target Names	PRM4
Protein Names	Recommended name: Pheromone-regulated membrane protein 4
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.