



Recombinant Escherichia coli Uncharacterized protein yfhR (yfhR)

Product Code	CSB-CF302872ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77538
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MALPVNKRVPKILFILFVVAFCVYLVPRVAINFFYYPDDKIYGPDPWSAESVEFT AKDGT RLQGWFISSSTGPADNAIATIIHAHG NAGNMSAHWPLVSWLPERNFNVFMFDY RGFGKSK GTPSQAGLLDDTQSAINVVRHRSDVNPQRLVLFQGQSIGGANILDVIGRGDREGI RAVILD STFASYATIANQMIPGSGYLLDESYSGENYIASVSPILLIHKADHVIPWQHS EKLYS LAKEPKRLILIPDGEHIDAFSDRHGDVYREQMVFILSALNPQN
Source	in vitro E.coli expression system
Target Names	yfhR
Protein Names	Recommended name: Uncharacterized protein yfhR
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.